


Product datasheet for **GA206981**

Mouse Tmem223 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem223 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 223
Symbol:	Tmem223
Components:	GA206981G1, Tmem223 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206981G2, Tmem223 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206981G3, Tmem223 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CQE2</u>
Synonyms:	0610006I08Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)