

## Product datasheet for **GA206998**

### Mouse Hmg20a activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | Hmg20a CRISPRa kit - CRISPR gene activation of mouse high mobility group 20A |
| Symbol:       | Hmg20a   |
| UniProt ID:   | <a href="#">Q9DC33</a>   |
| Synonyms:     | 1200004E06Rik; 5730490E10Rik; Hmgxb1; lbraf                                  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)      |

#### Product images:



[View online >](#)