

Product datasheet for **GA207005**

Mouse Rsrc1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Rsrc1 CRISPRa kit - CRISPR gene activation of mouse arginine/serine-rich coiled-coil 1
Symbol:	Rsrc1
Components:	GA207005G1, Rsrc1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207005G2, Rsrc1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207005G3, Rsrc1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9DBU6
Synonyms:	1200013F24Rik; SRrp53
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)