


Product datasheet for **GA207019**

Mouse Tmem107 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem107 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 107
Symbol:	Tmem107
Components:	GA207019G1, Tmem107 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207019G2, Tmem107 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207019G3, Tmem107 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CPV0</u>
Synonyms:	1110004B13Rik; 2810049P21Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)