

Product datasheet for **GA207058**

Mouse Slc25a23 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc25a23 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 25 (mitochondrial carrier; phosphate carrier), member 23
Symbol:	Slc25a23
Components:	GA207058G1, Slc25a23 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207058G2, Slc25a23 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207058G3, Slc25a23 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6GQS1
Synonyms:	2310067G05Rik; SCaMC-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)