

Product datasheet for **GA207074**

Mouse Psm5 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Psm5 CRISPRa kit - CRISPR gene activation of mouse proteasome (prosome, macropain) 26S subunit, non-ATPase, 5 |
| Symbol: | Psm5 |
| UniProt ID: | Q8BJY1 |
| Synonyms: | 1500032A03Rik; AW475925; mKIAA0072; S5b; W91691 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)