

Product datasheet for **GA207077**

Mouse Polr3k activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Polr3k CRISPRa kit - CRISPR gene activation of mouse polymerase (RNA) III (DNA directed) polypeptide K
Symbol:	Polr3k
UniProt ID:	Q9CQZ7
Synonyms:	12.3kDa; 1500004O14Rik; A1196849; AU014893; C11; RPC10; RPC11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)