

Product datasheet for **GA207080**

Mouse **Ttc23** activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ttc23 CRISPRa kit - CRISPR gene activation of mouse tetratricopeptide repeat domain 23
Symbol:	Ttc23
UniProt ID:	<u>Q8CHY7</u>
Synonyms:	1600012K10Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)