

Product datasheet for GA207091

Mouse Thap4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Thap4 CRISPRa kit - CRISPR gene activation of mouse THAP domain containing 4
Symbol:	Thap4
Components:	GA207091G1, Thap4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207091G2, Thap4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207091G3, Thap4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6P3Z3</u>
Synonyms:	2010320B01Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

