

Product datasheet for **GA207095**

Mouse Upf3a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Upf3a CRISPRa kit - CRISPR gene activation of mouse UPF3 regulator of nonsense transcripts homolog A (yeast)
Symbol:	Upf3a
Components:	GA207095G1, Upf3a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207095G2, Upf3a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207095G3, Upf3a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	2600001C03Rik; 4930546M19Rik; RENT3A; UPF3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)