


Product datasheet for **GA207108**

Mouse Pus3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pus3 CRISPRa kit - CRISPR gene activation of mouse pseudouridine synthase 3
Symbol:	Pus3
Components:	GA207108G1, Pus3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207108G2, Pus3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207108G3, Pus3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9J138
Synonyms:	2610020J05Rik; 5730412F04Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)