

Product datasheet for **GA207115**

Mouse Pgap4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem246 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 246
Symbol:	Pgap4
Components:	GA207115G1, Pgap4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207115G2, Pgap4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207115G3, Pgap4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q91YV9
Synonyms:	2810432L12Rik; 9330170P15Rik; AI835809
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)