

Product datasheet for **GA207117**

Mouse Polr3d activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Polr3d CRISPRa kit - CRISPR gene activation of mouse polymerase (RNA) III (DNA directed) polypeptide D
Symbol:	Polr3d
UniProt ID:	<u>Q91WD1</u>
Synonyms:	44kDa; 2810426M17Rik; AI326118; AW489084; BN51T; RPC4; TSNB51
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)