


Product datasheet for **GA207123**

Mouse Pi4k2b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pi4k2b CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol 4-kinase type 2 beta
Symbol:	Pi4k2b
Components:	GA207123G1, Pi4k2b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207123G2, Pi4k2b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207123G3, Pi4k2b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8CBQ5
Synonyms:	2610042N09Rik; 4933409G22Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)