

Product datasheet for GA207147

Mouse Dynlt3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dynlt3 CRISPRa kit - CRISPR gene activation of mouse dynein light chain Tctex-type 3
Symbol:	Dynlt3
Components:	GA207147G1, Dynlt3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207147G2, Dynlt3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207147G3, Dynlt3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P56387</u>
Synonyms:	2310075M16Rik; AU042796; Tctell
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

