


Product datasheet for **GA207174**

Mouse Sft2d3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sft2d3 CRISPRa kit - CRISPR gene activation of mouse SFT2 domain containing 3
Symbol:	Sft2d3
Components:	GA207174G1, Sft2d3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207174G2, Sft2d3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207174G3, Sft2d3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CSV6</u>
Synonyms:	2610206C24Rik; AI847965; AW046173
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)