

## Product datasheet for **GA207241**

### Mouse Washc3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Washc3 CRISPRa kit - CRISPR gene activation of mouse WASH complex subunit 3
Symbol:	Washc3
Components:	GA207241G1, Washc3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207241G2, Washc3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207241G3, Washc3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<a href="#">Q9CR27</a>
Synonyms:	2900091E11Rik; 5730495F03Rik; Ccdc53
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)