

Product datasheet for **GA207269**

Mouse Stk35 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Stk35 CRISPRa kit - CRISPR gene activation of mouse serine/threonine kinase 35
Symbol:	Stk35
UniProt ID:	<u>Q80ZW0</u>
Synonyms:	1700054C12Rik; CLIK1; Stk35I1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)