

Product datasheet for **GA207292**

Mouse 2210010C04Rik activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | 2210010C04Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 2210010C04 gene |
| Symbol: | 2210010C04Rik |
| Components: | GA207292G1, 2210010C04Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207292G2, 2210010C04Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207292G3, 2210010C04Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms: | AV072249 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)