

Product datasheet for **GA207309**

Mouse Txnrc8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Txnrc8 CRISPRa kit – CRISPR gene activation of mouse thioredoxin domain containing 8
Symbol:	Txnrc8
Components:	GA207309G1, Txnrc8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207309G2, Txnrc8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207309G3, Txnrc8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q69AB2</u>
Synonyms:	4930429J24Rik; Sptrx-3; Sptrx3; Trx6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

