

Product datasheet for **GA207311**

Mouse Soga3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Soga3 CRISPRa kit - CRISPR gene activation of mouse SOGA family member 3
Symbol:	Soga3
Components:	GA207311G1, Soga3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207311G2, Soga3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207311G3, Soga3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6NZL0</u>
Synonyms:	6330407J23Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)