

Product datasheet for GA207333

Mouse Dusp28 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dusp28 CRISPRa kit - CRISPR gene activation of mouse dual specificity phosphatase 28
Symbol:	Dusp28
Components:	GA207333G1, Dusp28 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA207333G2, Dusp28 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA207333G3, Dusp28 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q8BTR5</u>
Synonyms:	0710001B24Rik; AV005521
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

