

## Product datasheet for **GA207356**

### Mouse Eepd1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | Eepd1 CRISPRa kit - CRISPR gene activation of mouse endonuclease/exonuclease/phosphatase family domain containing 1                                                                                                                                                                                        |
| Symbol:       | Eepd1                                                                                                                                                                                                                                                                                                      |
| Components:   | GA207356G1, Eepd1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207356G2, Eepd1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207356G3, Eepd1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q3TGW2</a></u>                                                                                                                                                                                                                                                                              |
| Synonyms:     | 2310005P05Rik; AI481005                                                                                                                                                                                                                                                                                    |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)