

Product datasheet for **GA207357**

Mouse Dhx40 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dhx40 CRISPRa kit - CRISPR gene activation of mouse DEAH (Asp-Glu-Ala-His) box polypeptide 40
Symbol:	Dhx40
Components:	GA207357G1, Dhx40 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207357G2, Dhx40 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207357G3, Dhx40 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PE54</u>
Synonyms:	2410016C14Rik; ARG147; DDX40; PAD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)