

Product datasheet for **GA207362**

Mouse Tmem167b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem167b CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 167B
Symbol:	Tmem167b
Components:	GA207362G1, Tmem167b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207362G2, Tmem167b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207362G3, Tmem167b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80X45</u>
Synonyms:	2010200O16Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)