


## Product datasheet for **GA207370**

### Mouse Tvp23b activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Tvp23b CRISPRa kit - CRISPR gene activation of mouse trans-golgi network vesicle protein 23B                                                                                                                                                                                                                  |
| Symbol:       | Tvp23b                                                                                                                                                                                                                                                                                                        |
| Components:   | GA207370G1, Tvp23b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207370G2, Tvp23b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207370G3, Tvp23b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q9D8T4</a></u>                                                                                                                                                                                                                                                                                 |
| Synonyms:     | 1810036I24Rik; Fam18b; Fam18b1                                                                                                                                                                                                                                                                                |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)