

Product datasheet for **GA207393**

Mouse Pigm activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pigm CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol glycan anchor biosynthesis, class M
Symbol:	Pigm
UniProt ID:	Q8C2R7
Synonyms:	4933437L05Rik; AI644422; AV302080; C920011G20Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)