


Product datasheet for **GA207396**

Mouse Mgat4c activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mgat4c CRISPRa kit - CRISPR gene activation of mouse MGAT4 family, member C
Symbol:	Mgat4c
Components:	GA207396G1, Mgat4c gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207396G2, Mgat4c gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207396G3, Mgat4c gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9D306
Synonyms:	9130411I17Rik; GntIVh
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)