

Product datasheet for GA207437

Mouse 4930522H14Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	4930522H14Rik CRISPRa kit – CRISPR gene activation of mouse RIKEN cDNA 4930522H14 gene
Symbol:	4930522H14Rik
Components:	GA207437G1, 4930522H14Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207437G2, 4930522H14Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207437G3, 4930522H14Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9CPZ3
Synonyms:	1700007J23Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

