

Product datasheet for **GA207468**

Mouse Reg4 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Reg4 CRISPRa kit - CRISPR gene activation of mouse regenerating islet-derived family, member 4 |
| Symbol: | Reg4 |
| UniProt ID: | Q9D8G5 |
| Synonyms: | 2010002L15Rik; GISP; RELP |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)