

Product datasheet for **GA207471**

Mouse Slc25a37 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc25a37 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 25, member 37
Symbol:	Slc25a37
Components:	GA207471G1, Slc25a37 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207471G2, Slc25a37 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207471G3, Slc25a37 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q920G8
Synonyms:	1700020E22Rik; 4930513O14Rik; 4930526G11Rik; AI848481; C330015G08Rik; frascati; Mfrn
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)