

Product datasheet for **GA207475**

Mouse 2310057J18Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	2310057J18Rik CRISPRa kit – CRISPR gene activation of mouse RIKEN cDNA 2310057J18 gene
Symbol:	2310057J18Rik
Components:	GA207475G1, 2310057J18Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207475G2, 2310057J18Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207475G3, 2310057J18Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8C6C9</u>
Synonyms:	2210414P04Rik; Leg1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

