

Product datasheet for **GA207490**

Mouse Slc48a1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc48a1 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 48 (heme transporter), member 1
Symbol:	Slc48a1
Components:	GA207490G1, Slc48a1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207490G2, Slc48a1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207490G3, Slc48a1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9D8M3
Synonyms:	4930570C03Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)