


Product datasheet for **GA207498**

Mouse Ddx47 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ddx47 CRISPRa kit - CRISPR gene activation of mouse DEAD (Asp-Glu-Ala-Asp) box polypeptide 47
Symbol:	Ddx47
Components:	GA207498G1, Ddx47 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207498G2, Ddx47 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207498G3, Ddx47 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CWX9</u>
Synonyms:	4930588A18Rik; C77285
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)