

Product datasheet for **GA207518**

Mouse Dalrd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dalrd3 CRISPRa kit - CRISPR gene activation of mouse DALR anticodon binding domain containing 3
Symbol:	Dalrd3
Components:	GA207518G1, Dalrd3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207518G2, Dalrd3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207518G3, Dalrd3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PJN8</u>
Synonyms:	6330580J24Rik; C77829; mir-425
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)