

Product datasheet for **GA207531**

Mouse Snap47 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Snap47 CRISPRa kit - CRISPR gene activation of mouse synaptosomal-associated protein, 47
Symbol:	Snap47
UniProt ID:	Q8R570
Synonyms:	1110031B06Rik; SNAP-47
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)