

Product datasheet for GA207535

Mouse Wdr83 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr83 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain containing 83
Symbol:	Wdr83
Components:	GA207535G1, Wdr83 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207535G2, Wdr83 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207535G3, Wdr83 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9DAJ4
Synonyms:	1500041N16Rik; Morg1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)