


Product datasheet for **GA207544**

Mouse Ddx55 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ddx55 CRISPRa kit - CRISPR gene activation of mouse DEAD (Asp-Glu-Ala-Asp) box polypeptide 55
Symbol:	Ddx55
Components:	GA207544G1, Ddx55 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207544G2, Ddx55 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207544G3, Ddx55 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZPL9</u>
Synonyms:	2810021H22Rik; mKIAA1595
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)