

## Product datasheet for **GA207550**

### Mouse **S100a16** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	S100a16 CRISPRa kit - CRISPR gene activation of mouse S100 calcium binding protein A16
Symbol:	S100a16
UniProt ID:	<a href="#">Q9D708</a>
Synonyms:	2300002L21Rik; AI325039; AI663996; DT1P1A7; S100F
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)