

Product datasheet for GA207570

Mouse Saraf activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Saraf CRISPRa kit – CRISPR gene activation of mouse store-operated calcium entry-associated regulatory factor
Symbol:	Saraf
Components:	GA207570G1, Saraf gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207570G2, Saraf gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207570G3, Saraf gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8R3Q0</u>
Synonyms:	1810045K07Rik; Tmem66
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

