

Product datasheet for **GA207584**

Mouse Plpp3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plpp3 CRISPRa kit - CRISPR gene activation of mouse phospholipid phosphatase 3
Symbol:	Plpp3
UniProt ID:	Q99JY8
Synonyms:	1110003O22Rik; 2610002D05Rik; AV025606; D4Bwg0538e; D4Bwg1535e; Lpp3; Ppab2b; Ppap2b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)