


## Product datasheet for **GA207615**

### Mouse Ooep activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | Ooep CRISPRa kit - CRISPR gene activation of mouse oocyte expressed protein                                                                                                                                                                                                                             |
| Symbol:       | Ooep                                                                                                                                                                                                                                                                                                    |
| Components:   | GA207615G1, Ooep gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207615G2, Ooep gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207615G3, Ooep gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q9CWE6</a>                                                                                                                                                                                                                                                                                  |
| Synonyms:     | 2410146L05Rik; Floped; Moep19; Sddr                                                                                                                                                                                                                                                                     |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)