


Product datasheet for **GA207616**

Mouse Tppp3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tppp3 CRISPRa kit - CRISPR gene activation of mouse tubulin polymerization-promoting protein family member 3
Symbol:	Tppp3
Components:	GA207616G1, Tppp3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207616G2, Tppp3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207616G3, Tppp3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CRB6</u>
Synonyms:	2700055K07Rik; Ceacam9; CGI-38; mmCGM8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)