

Product datasheet for **GA207631**

Mouse Retreg3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Retreg3 CRISPRa kit - CRISPR gene activation of mouse reticulophagy regulator family member 3
Symbol:	Retreg3
Components:	GA207631G1, Retreg3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207631G2, Retreg3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207631G3, Retreg3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9CQV4
Synonyms:	1300010M03Rik; 4933404C01Rik; AI551748
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)