

## Product datasheet for **GA207640**

### Mouse Pdf activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                      |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                         |
| Description:  | Pdf CRISPRa kit - CRISPR gene activation of mouse peptide deformylase (mitochondrial)                                                                                                                                                                                                                |
| Symbol:       | Pdf                                                                                                                                                                                                                                                                                                  |
| Components:   | GA207640G1, Pdf gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207640G2, Pdf gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207640G3, Pdf gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms:     | 2610019N19Rik                                                                                                                                                                                                                                                                                        |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                              |

#### Product images:



[View online »](#)