


Product datasheet for **GA207643**

Mouse Tmem178 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem178 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 178
Symbol:	Tmem178
Components:	GA207643G1, Tmem178 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207643G2, Tmem178 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207643G3, Tmem178 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CZ16</u>
Synonyms:	2810417M05Rik; Tmem178a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)