

Product datasheet for **GA207661**

Mouse Dmac2l activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dmac2l CRISPRa kit - CRISPR gene activation of mouse distal membrane arm assembly complex 2 like
Symbol:	Dmac2l
Components:	GA207661G1, Dmac2l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207661G2, Dmac2l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207661G3, Dmac2l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CRA7</u>
Synonyms:	1110015E18Rik; ATPW
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)