

Product datasheet for **GA207664**

Mouse Mrnip activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mrnip CRISPRa kit - CRISPR gene activation of mouse MRN complex interacting protein
Symbol:	Mrnip
Components:	GA207664G1, Mrnip gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207664G2, Mrnip gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207664G3, Mrnip gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D1F5</u>
Synonyms:	2900016G09Rik; 3010026O09Rik; AI851514; AW475975
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)