

Product datasheet for GA207671

Mouse Gpn3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Gpn3 CRISPRa kit – CRISPR gene activation of mouse GPN-loop GTPase 3
Symbol:	Gpn3
Components:	GA207671G1, Gpn3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA207671G2, Gpn3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA207671G3, Gpn3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q9D3W4</u>
Synonyms:	A930018B01Rik; Atpbd1c; D5Ert708e
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

