

## Product datasheet for **GA207673**

### Mouse *Arc4* activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | <i>Arc4</i> CRISPRa kit - CRISPR gene activation of mouse actin related protein 2/3 complex, subunit 4 |
| Symbol:       | <i>Arc4</i>  |
| UniProt ID:   | <a href="#">P59999</a>   |
| Synonyms:     | 20kDa; 5330419I20Rik; AI327076; p20-Arc  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                |

#### Product images:



[View online >](#)